



Regulator, InGage



DR8033-7003P4



DR8133-7003P4



DR8333-7003P4

The Flotec **InGage™ Regulator** is a medical gas regulator with a unique internal pressure gauge built into the body of the regulator with a window to show the cylinder pressure. The gauge features a color coded 0-3000 PSI range that can be reset in the field. The all brass, high pressure, yoke connector has a 10-micron sintered brass inlet filter and a permanent O-ring type yoke washer seal. The regulator is a single stage, piston type, backpressure compensated, pressure reducer with an integral relief valve and a “click style” flowmeter that has 12 flow positions that will operate in any orientation or during extreme vibration. The unique flow rotor assembly has a second integrated 20-micron, sintered stainless steel filter built into the laser drilled orifice plate providing an absolute level of contamination protection. Currently available in nine different Medical Gases, 16 colors, six pressure settings, five International standards, two different metals and a nearly unlimited number of outlet connections. The **InGage™ Regulator** design allows the user to customize the regulator so that it will fit almost any requirement. The unique and distinctive Flotec Dodecagon (12 sided) shape helps prevent confusion with other brands, prevents the regulator from rolling off tabletops, and makes it easier to crack pecans, kill bugs and drive nails.

- Unique Internal Pressure Gauge
- Unique Field Adjustable Pressure Gauge
- Unique Laser Drilled Orifice Plate with Integral Filter
- Unique Dome Handle
- Unique O-Ring Type Yoke Washer Seal
- Unique Swivel Yoke & Body
- Unique Dodecagon (12 sided) shape
- Unique Double Filter System
- Unique Thread Insert Design
- Unique Custom Laser Etching available
- MRI Conditional
- 9 Medical Gas configurations available in 13 Different Anodized Colors
- All Metal Construction (available in Anodized Aluminum and Plated Brass)
- FDA CGMP (CFR 21, Parts 803-820), ISO9001:2008, and ISO13485:2003
- **Manufactured exclusively in the United States of America**



375 PSI

2000 PSI

Yoke Washer Seal

Dome Handle



Part Number	Description
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb
DR8133-7003P4	CGA 870, One DISS 1240 Check Valve, 0-25 LPM, Hose Barb
DR8333-7003P4	CGA 870, Two DISS 1240 Check Valves, 0-25 LPM, Hose Barb
DR5103-6003	CGA 540, One DISS Check Valve, Pressure Reducer, Plain Cap
DR5303-6003	CGA 540, Two DISS Check Valves, Pressure Reducer, Plain Cap
DR5303-7003P4	CGA 540, Two DISS 1240 Check Valves, 0-25 LPM, Hose Barb

Match Ambulance outlets with Regulator outlets for quicker patient transfers

- Ohio Quick Disconnect and DISS-1240 Check Valve
- Chemetron Quick Disconnect and DISS-1240 Check Valve
- DISS-1240 Check Valve



Part Number	Description
DR8E33-7003P4	CGA 870, One Ohio & One Check Valve, 0-25 LPM, Hose Barb
DR8Q33-7003P4	CGA 870, One Chemetron & One Check Valve, 0-25 LPM, Hose Barb

Anodized Aluminum

Part Number	Description	
DR8032-7002P4	CGA 870, 0-25 LPM, Hose Barb, Silver (Clear)	
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb, Green	
DR8034-7004P4	CGA 870, 0-25 LPM, Hose Barb, Yellow	
DR8035-7005P4	CGA 870, 0-25 LPM, Hose Barb, Blue	
DR8036-7006P4	CGA 870, 0-25 LPM, Hose Barb, Red	
DR8037-7007P4	CGA 870, 0-25 LPM, Hose Barb, Brown	
DR8038-7008P4	CGA 870, 0-25 LPM, Hose Barb, Gold	
DR8039-7009P4	CGA 870, 0-25 LPM, Hose Barb, Black	
DR803A-700AP4	CGA 870, 0-25 LPM, Hose Barb, Purple	
DR803B-700BP4	CGA 870, 0-25 LPM, Hose Barb, Orange	
DR803C-700CP4	CGA 870, 0-25 LPM, Hose Barb, Camouflage Desert	
DR803P-700PP4	CGA 870, 0-25 LPM, Hose Barb, Pink	
DR803W-700WP4	CGA 870, 0-25 LPM, Hose Barb, Camouflage Woodland	

Plated Brass

Part Number	Description
DR803T-700TP4	CGA 870, 0-25 LPM, Hose Barb, Silver (Tri-Metal)
DR8033-7003P4	CGA 870, 0-25 LPM, Hose Barb, 23k Gold Plated

