

----- FLOW RANGE:

See B-Flows chart for part number and list price

----- COLOR OPTION:

- NONE - Lt Blue Anodized Body
- 2 - Silver Anodized Body
- 3 - Green Anodized Body
- 4 - Yellow Anodized Body
- 5 - Blue Anodized Body
- 6 - Red Anodized Body
- 7 - Brown Anodized Body
- 8 - Gold Anodized Body
- 9 - Black Anodized Body
- A - Purple Anodized Body
- B - Orange Anodized Body
- G - Gold Plated Brass Body (No Aluminum)
- T - Tri-Metal Plated Brass Body (No Aluminum)

-OR-

FITTING OUTLET:

- A - 1 LPM Orifice
- B - 2 LPM Orifice
- C - 3 LPM Orifice
- D - 4 LPM Orifice
- E - 5 LPM Orifice
- F - 6 LPM Orifice
- H - 8 LPM Orifice
- J - 10 LPM Orifice
- L - 12 LPM Orifice
- M - 15 LPM Orifice
- N - 20 LPM Orifice
- P - 25 LPM Orifice
- T - 60 LPM Orifice
- V - 70 LPM Orifice
- W - 1/2 LPM Orifice (discontinued)

----- FLOW CARTRIDGE OPTIONS:

- 0 - No Flow Cartridge Options

-OR-

CAP OPTIONS:

- 0 - Plain Cap, No Outlets
- 1 - Hose Barb 1/8" MNPT
- 2 - 1/8" FNPT
- 4 - 3/16" Female SAE Boss Port
- 6 - Hose Barb 3/16" SAE Boss
- 7 - DISS-1240 Fitting
- B - One DISS-1240 Check Valve
- H - Hansen QD
- P - Puritain Bennett QD

----- WINDOW & CONTROLLED FLOW OUTLET POSITIONS:

- 0 - Window @ 0 & Hose Barb @ 180
- 1 - Window @ 0 & DISS-1240 Fitting @ 180
- 2 - Window @ 0 & 1/8" Female Pipe Port @ 180
- 7 - Window @ 0 & Hose Barb @ 0
- H - Window @ 0 & Clear Hose Barb @ 180
- N - Window @ 0 & DISS-1240 Fitting & HBD2 @ 180
- W - One DISS-1180 Check Valve

-OR-

CAP OUTLET POSITIONS:

- 0 - Outlet on Centerline
- 2 - Outlet @ 0
- 3 - Outlet @ 90
- 6 - Outlet @ 180
- 9 - Outlet @ 270

----- FLOW CARTRIDGE CONFIGURATION:

- 6 - Cap
- 7 - Large Diameter Knurled Knob
- 8 - Shrouded Knob

----- COLOR OPTION:

- NONE - Lt Blue Anodized Body
- 2 - Silver Anodized Body
- 3 - Green Anodized Body
- 4 - Yellow Anodized Body
- 5 - Blue Anodized Body
- 6 - Red Anodized Body
- 7 - Brown Anodized Body
- 8 - Gold Anodized Body
- 9 - Black Anodized Body
- A - Purple Anodized Body
- B - Orange Anodized Body
- G - Gold Plated Brass Body (No Aluminum)
- T - Tri-Metal Plated Brass Body (No Aluminum)

----- BODY OPTION(S):

- 0 - No Body Options (\$0)

----- PRESSURE REDUCED BODY OUTLET(S):

- 0 - No Pressure Body Outlets
- 1 - One DISS-1240 Check Valve
- 2 - One 1/8" FNPT (Plugged)
- 3 - Two DISS-1240 Check Valves
- 4 - Two 1/8" FNPT (Plugged)
- 6 - One DISS-1240 Check Valve & One 1/8" FNPT (Plugged)
- 7 - One DISS-1160 Check Valve @ 180
- B - One DISS-1240 Check Valve @ 180
- E - One Female Ohio QD & One DISS-1240 Check Valve
- H - One Female Ohio QD
- J - One DISS-1240 Check Valve @ 330
- P - Two DISS-1240 Check Valves w/ One HBD2-50-10
- Q - One Female Chemetron QD & One DISS-1240 Check Valve
- R - Hansen QD
- S - One DISS-1060 Check Valve @ 180
- U - One Hansen QD & One DISS-1240 Check Valve
- V - One 1/8" FNPT (Plugged) @ 180
- W - One DISS-1180 Check Valve

----- CYLINDER SIZE & TYPE:

- XX - No Cylinder
- CB - B size Aluminum
- CC - C size Aluminum
- CD - D size Aluminum
- CJ - Jumbo D size Aluminum
- CE - E size Aluminum
- CM - M size Aluminum
- CFA - A size Composite
- CFB - B size Composite
- CFC - C size Composite
- CFD - D size Composite
- CHB - B size 3000 PSI Aluminum
- CHC - C size 3000 PSI Aluminum
- CHD - D size 3000 PSI Aluminum
- CHE - E size 3000 PSI Aluminum
- CLD - D size Lightweight Steel
- CLE - E size Lightweight Steel

----- VALVE OPTIONS:

- 1 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Welded (replaced by T)
- 3 - Integral Shut Off Valve w/ Knurled Knob, Bent Handle, Adapter (replaced by B)
- 4 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Adapter (discontinued)
- 5 - Integral Shut Off Valve w/ Knurled Knob, Bent Handle, CGA 280, Adapter (replaced by C)
- 6 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, CGA 346, Welded (replaced by D)
- 7 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, 1/2-14 NGT, Adapter (discontinued)
- 8 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, 1/2-14 NGT, Welded (replaced by G)
- 9 - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, External 1.4", 3000 PSI Gauge, Welded
- A - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, 1/2-14 NGT, External 1.4", 3000 PSI Gauge, Welded (discontinued)
- P - Integral Shut Off Valve w/ Knurled Knob, No Handle (replaced by H)

OutGage

- B - Integral Shut Off Valve w/ Knurled Knob, Bent Handle, Adapter, External 1.25", 3000 PSI Gauge
- C - Integral Shut Off Valve w/ Knurled Knob, Bent Handle, CGA 280, Adapter, External 1.25", 3000 PSI Gauge
- D - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, CGA 346, Welded, External 1.25", 3000 PSI Gauge
- G - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, 1/2-14 NGT, Welded, External 1.25", 3000 PSI Gauge
- H - Integral Shut Off Valve w/ Knurled Knob, No Handle, External 1.25", 3000 PSI Gauge
- T - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Welded, External 1.25", 3000 PSI Gauge

RECONDITIONED

- E - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Welded, External 1.25" 4000 PSI Gauge, for 3000 psi cylinders
- F - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Adapter, External 1.25" 4000 PSI Gauge, for 3000 psi cylinders
- J - Integral Shut Off Valve w/ Knurled Knob, No Handle, External 1.25", 3000 PSI Gauge
- K - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Welded, External 1.25" 3000 PSI Gauge
- L - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Adapter, External 1.25" 3000 PSI Gauge
- R - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Welded
- S - Integral Shut Off Valve w/ Knurled Knob, Halo Handle, Adapter

----- PRODUCT GROUP:

FP - Flopac cylinder mounted Ingage Regulator