

---- Fitting end

- 12 - 1/8" MNPT
- 14 - 1/4" MNPT
- 22 - 1/8" FNPT
- 2S - 1/8" FNPT Swivel
- DCV - DISS Check Valve
- DF - DISS Fitting
- HB-4 - 1/4" Hose Barb
- HB-5 - 5/16" Hose Barb
- HN - Hex Nut Female
- KN - Knurled Nut Female
- 90-HB-4 - 90 Elbow 1/4" Hose Barb

---- QD End

- Quick Disconnect Female Chemetron Carbon Dioxide - U2
- Quick Disconnect Female Chemetron Medical Air - D2
- Quick Disconnect Female Chemetron Nitrous Oxide - N2
- Quick Disconnect Female Chemetron Oxygen - C2
- Quick Disconnect Female Chemetron Vacuum (Suction) - V2
- Quick Disconnect Female Chemetron WAGD (EVAC) - E2
- Quick Disconnect Female Ohio Carbon Dioxide - Y2
- Quick Disconnect Female Ohio Medical Air - A2
- Quick Disconnect Female Ohio Nitrous Oxide - H2
- Quick Disconnect Female Ohio Oxygen - O2
- Quick Disconnect Female Ohio Vacuum (Suction) - J2
- Quick Disconnect Female Puritan Bennet Nitrous Oxide - K2
- Quick Disconnect Female Puritan Bennet Vacuum (Suction) - M2
- Quick Disconnect Female Puritan Bennet WAGD (EVAC) - B2
- Quick Disconnect Female Puritan Bennett Medical Air - R2
- Quick Disconnect Female Puritan Bennett Oxygen - P2
- Quick Disconnect Female Schrader Medical Air - T2
- Quick Disconnect Female Schrader Nitrogen - Q2
- Quick Disconnect Female Schrader Nitrous Oxide - G2
- Quick Disconnect Female Schrader Oxygen - S2
- Quick Disconnect Female Schrader Vacuum (Suction) - W2
- Quick Disconnect Hex Nut Female (NIST) National Institute Standards & Technology Medical Air - D8
- Quick Disconnect Hex Nut Female (NIST) National Institute Standards & Technology Oxygen - D6
- Quick Disconnect Knurled Nut Female 90 degree elbow (AFNOR) French Medical Air - H8
- Quick Disconnect Knurled Nut Female 90 degree elbow (AFNOR) French Oxygen - H6
- Quick Disconnect Knurled Nut Female (SIS) Australian Medical Air - K8
- Quick Disconnect Knurled Nut Female (SIS) Australian Oxygen - K6
- Quick Disconnect Male Chemetron Carbon Dioxide - U1
- Quick Disconnect Male Chemetron Medical Air - D1
- Quick Disconnect Male Chemetron Nitrous Oxide - N1
- Quick Disconnect Male Chemetron Oxygen - C1
- Quick Disconnect Male Chemetron Vacuum (Suction) - V1
- Quick Disconnect Male Chemetron WAGD (EVAC) - E1
- Quick Disconnect Male Medstar Medical Air x 1/8" MNPT - L3
- Quick Disconnect Male Medstar Nitrous Oxide - L1
- Quick Disconnect Male Medstar Oxygen x 1/8" MNPT - L5
- Quick Disconnect Male Medstar Vacuum (Suction) x 1/8" MNPT - L7
- Quick Disconnect Male non-swivel Schrader Medical Air - T1
- Quick Disconnect Male non-swivel Schrader Nitrogen - Q1
- Quick Disconnect Male non-swivel Schrader Nitrous Oxide - G1
- Quick Disconnect Male non-swivel Schrader Oxygen - S1
- Quick Disconnect Male non-swivel Schrader Vacuum (Suction) - W1
- Quick Disconnect Male Ohio Carbon Dioxide - Y1
- Quick Disconnect Male Ohio Medical Air - A1
- Quick Disconnect Male Ohio Medical Air (ISO Black/White) - A1BW
- Quick Disconnect Male Ohio Nitrous Oxide - H1
- Quick Disconnect Male Ohio Oxygen - O1
- Quick Disconnect Male Ohio Oxygen (ISO White) - O1W
- Quick Disconnect Male Ohio Vacuum (Suction) - J1
- Quick Disconnect Male Ohio Vacuum (Suction) (ISO Yellow) - J1Y
- Quick Disconnect Male Ohio WAGD (EVAC) - F1
- Quick Disconnect Male Ohio WAGD (EVAC) - F2
- Quick Disconnect Male Oxequip Medical Air - Z3
- Quick Disconnect Male Oxequip Nitrous Oxide - Z1
- Quick Disconnect Male Oxequip Oxygen - Z5
- Quick Disconnect Male Oxequip Vacuum (Suction) - Z7
- Quick Disconnect Male Puritan Bennet Nitrous Oxide K1
- Quick Disconnect Male Puritan Bennet Vacuum (Suction) - M1
- Quick Disconnect Male Puritan Bennet Vacuum (Suction) (ISO Yellow) - M1Y
- Quick Disconnect Male Puritan Bennet WAGD (EVAC) - B1
- Quick Disconnect Male Puritan Bennett Medical Air - R1
- Quick Disconnect Male Puritan Bennett Medical Air (ISO Yellow) - R1Y
- Quick Disconnect Male Puritan Bennett Oxygen - P1
- Quick Disconnect Male Puritan Bennett Oxygen (ISO White) - P1W
- Quick Disconnect Male swivel Schrader Medical Air - T3
- Quick Disconnect Male swivel Schrader Nitrogen - Q3
- Quick Disconnect Male swivel Schrader Nitrous Oxide - G3
- Quick Disconnect Male swivel Schrader Oxygen - S3
- Quick Disconnect Male swivel Schrader Vacuum (Suction) - W3
- Quick Disconnect Male (BS) British Standard Medical Air - B8
- Quick Disconnect Male (BS) British Standard Oxygen - B6
- Quick Disconnect Male (DIN) Germany Medical Air - F8
- Quick Disconnect Male (DIN) Germany Oxygen - F6

----- PRODUCT GROUP:
 QD - Quick Disconnect Fitting