



## Medical Oxygen Cylinders

### Aluminum

2015 or 2216 psi Rating (139 bar)

Manufactured to DOT and Transport Canada standards. Non-magnetic.  
0.750-16 UNF thread size up to Size CM which is 1.125-12 UNF.  
Sizes B, C, D, Jumbo D and E



### Composite Carbon-Fiber

3000 psi Rating (207 bar)

Wrapped in SE Fiberglass over Aluminum shell and finished with either Green or Grey tint and clear resin. 0.750-16 UNF Thread Size.  
Sizes CFA, CFB, CFC and CFD

### High Pressure Aluminum

3000 psi Rating (207 bar)


Higher-Strength aluminum alloy with minimal weight increase  
Sizes B, C, D and E

**FloGripp!**

**New Post Valve Replacement Knob**

Available for all CGA-870 Post Valve applications

**See reverse side for all available choices!**





## Composite (Carbon-Fiber) Oxygen Cylinders

### Size FB composite cylinders (3.40" diameter)

244 volume in liters (8.6 cu ft) 3000 psi (207 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CFB520BH	CGA 540 vertical	None	2.3 lbs.	13.95"
CFB521BH	CGA 540 vertical	RH	2.4 lbs.	13.95"
CFB800BH	CGA 870	None	2.4 lbs.	15.58"
CFB802BH	CGA 870	LH	2.4 lbs.	15.58"
CFB810BH	CGA 870 w/ toggle	None	2.6 lbs.	15.70"
CFB830BH	CGA-870 FloGripp	None	2.6 lbs.	16.60"

### Size FC composite cylinders (4.00" diameter)

437 volume in liters (15.4 cu ft) 3000 psi (207 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CFC520BH	CGA 540 vertical	None	3.3 lbs.	16.45"
CFC521BH	CGA 540 vertical	RH	3.4 lbs.	16.45"
CFC800BH	CGA 870	None	3.4 lbs.	18.08"
CFC802BH	CGA 870	LH	3.4 lbs.	18.08"
CFC810BH	CGA 870 w/ toggle	None	3.6 lbs.	18.20"
CFC830BH	CGA-870 FloGripp	None	3.6 lbs.	19.23"

### Size FD composite cylinders (4.00" diameter)

613 volume in liters (21.6 cu ft) 3000 psi (207 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CFD520BH	CGA 540 vertical	None	4.5 lbs.	17.25"
CFD521BH	CGA 540 vertical	RH	4.6 lbs.	17.25"
CFD800BH	CGA 870	None	4.6 lbs.	18.88"
CFD802BH	CGA 870	LH	4.6 lbs.	18.88"
CFD810BH	CGA 870 w/ toggle	None	4.8 lbs.	19.00"
CFD830BH	CGA-870 FloGripp	None	4.8 lbs.	20.03"

## Aluminum Oxygen Cylinders

### Size B aluminum cylinders (3.20" diameter)

170 volume in liters (6.0 cu ft) 2216 psi (153 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CB520B	CGA 540 vertical	None	2.9lbs.	13.55"
CB521B	CGA 540 vertical	RH	2.9lbs.	13.55"
CB800B	CGA 870	None	2.9 lbs.	15.18"
CB802B	CGA 870	LH	2.9 lbs.	15.18"
CB810B	CGA 870 w/ toggle	None	3.0 lbs.	15.30"
CB830B	CGA-870 FloGripp	None	3.0 lbs.	16.30"

### Size C aluminum cylinders (4.37" Diameter)

255 volume in liters (9 cu ft) 2015 psi (139 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CC520B	CGA 540 vertical	None	4.3lbs.	12.75"
CC521B	CGA 540 vertical	RH	4.3lbs.	12.75"
CC800B	CGA 870	None	4.3 lbs.	14.32"
CC802B	CGA 870	LH	4.3 lbs.	14.32"
CC810B	CGA 870 w/ toggle	None	4.4 lbs.	14.41"
CC830B	CGA-870 FloGripp	None	4.4 lbs.	15.44"

### Size D aluminum cylinders (4.37" diameter)

425 volume in liters (15 cu ft) 2015 psi (139 bar)

Part Number	CGA-Valve	Gauge	Weight	Height
CD520B	CGA 540 vertical	None	5.8 lbs.	18.25"
CD521B	CGA 540 vertical	RH	5.8 lbs.	18.25"
CD800B	CGA 870	None	6.2 lbs.	19.88"
CD802B	CGA 870	LH	6.2 lbs.	19.88"
CD810B	CGA 870 w/ toggle	None	6.4 lbs.	20.00"
CD812B	CGA 870 w/ toggle	LH	6.4 lbs.	20.0"
CD832B	CGA-870 FloGripp	None	6.4 lbs.	21.06"



Size Jumbo D aluminum cylinders (5.30" diameter)					
		650 volume in liters (22.6 cu ft)	2216 psi (153 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CJ520B	CGA 540 vertical		None	9.6 lbs.	18.88"
CJ521B	CGA 540 vertical		RH	9.7 lbs.	18.88"
CJ800B	CGA 870		None	9.7 lbs.	20.44"
CJ802B	CGA 870		LH	9.7 lbs.	20.44"
CJ810B	CGA 870 w/ toggle		None	9.9 lbs.	20.53"
CJ830B	CGA-870 FloGripp		None	9.9 lbs.	21.56"

Aluminum Oxygen Cylinders (Cont.)					
Size E aluminum cylinders (4.37" diameter)					
		680 volume in liters (24.2 cu ft)	2015 psi (139 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CE520B	CGA 540 vertical		None	8.0 lbs.	26.75"
CE521B	CGA 540 vertical		RH	8.1 lbs.	26.75"
CE800B	CGA 870		None	8.1 lbs.	28.38"
CE802B	CGA 870		LH	8.2 lbs.	28.38"
CE810B	CGA 870 w/ toggle		None	8.2 lbs.	28.50"
CE830B	CGA-870 FloGripp		None	8.2 lbs.	30.20"
Size M aluminum cylinders (8.00" diameter)					
		3469 volume in liters (122 cu ft)	2216 psi (139 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CM520B	CGA 540 vertical		None	40.8 lbs.	36.97"
CM521B	CGA 540 vertical		RH	40.9 lbs.	36.97"
CM800B	CGA 870		None	40.9 lbs.	5.88"
CM802B	CGA 870		LH	40.9 lbs.	35.88"
CM810B	CGA 870 w/ toggle		None	41.1 lbs.	36.00"
Aluminum High Pressure Oxygen Cylinders					
Size B HP aluminum cylinders (3.20" diameter)					
		233 volume in liters (8.2 cu ft)	3000 psi (207 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CHB520B	CGA 540 vertical		None	2.9 lbs.	13.55"
CHB521B	CGA 540 vertical		RH	2.9 lbs.	13.55"
CHB800B	CGA 870		None	2.9 lbs.	15.18"
CHB802B	CGA 870		LH	2.9 lbs.	15.18"
CHB810B	CGA 870 w/ toggle		None	3.0 lbs.	15.30"
CHB830B	CGA-870 FloGripp		None	3.0 lbs.	16.33"
Size C HP aluminum cylinders (4.40" diameter)					
		354 volume in liters (12.5 cu ft)	3000 psi (207 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CHC520B	CGA 540 vertical		None	4.7 lbs.	12.35"
CHC521B	CGA 540 vertical		RH	4.7 lbs.	12.35"
CHC800B	CGA 870		None	4.7 lbs.	13.98"
CHC802B	CGA 870		LH	4.7 lbs.	13.98"
CHC810B	CGA 870 w/ toggle		None	4.8 lbs.	14.10"
CHC830B	CGA-870 FloGripp		None	4.8 lbs.	15.13"
Size D HP aluminum cylinder (4.40" diameter)					
		605 volume in liters (21.3 cu ft)	3000 psi (207 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CHD520B	CGA 540 vertical		None	6.4 lbs.	18.05"
CHD521B	CGA 540 vertical		RH	6.4 lbs.	18.05"
CHD800B	CGA 870		None	6.4 lbs.	19.68"
CHD802B	CGA 870		LH	6.4 lbs.	19.68"
CHD810B	CGA 870 w/ toggle		None	6.5 lbs.	19.80"
CHD812B	CGA 870 w/ toggle		LH	6.5 lbs.	19.80"
CHD832B	CGA-870 FloGripp		None	6.5 lbs.	20.93"

Size E HP aluminum cylinder (4.40" diameter)					
		992 volume in liters (35.0 cu ft)	3000 psi (207 bar)		
Part Number	CGA-Valve		Gauge	Weight	Height
CHE520B	CGA 540 vertical		None	9.0 lbs.	27.05"
CHE521B	CGA 540 vertical		RH	9.0 lbs.	27.05"
CHE800B	CGA 870		None	9.0 lbs.	28.68"
CHE802B	CGA 870		LH	9.0 lbs.	28.68"
CHE810B	CGA 870 w/ toggle		None	9.1 lbs.	28.80"
CHE830B	CGA-870 FloGripp		None	9.1 lbs.	30.83"
Steel High Pressure Cylinders					
Size H Steel Cylinders (9.25" diameter)					
Part Number	CGA-Valve	Capacity		Weight	Height
CSH540B	CGA-540	2216 psi, 250 cf		110	51"
CSH577B	CGA-577	3600 psi, 387 cf		144	51"
CSH577BH	CGA-577	4500 psi, 444 cf		144	51"
CSH701B	CGA-701	6000 psi, 509 cf		188	51"