



Why "Single" when you can "T" or "Y"!

# Double Flow Control Valve Assemblies

## Y-Blocks



**All-Metal-Construction  
 No Plastic!**



Dual -Flow Control Valves from a single wall outlet.


Optional Pediatric and High Flow Control Valves

Optional removable Resident Nipple Adapter for CPAP of Humidifier applications

Optional DISS-1240 Check Valve, Female Ohio or Chemetron

Optional T or Y styles.


Available flows from 1/32nd of a liter thru 60 LPM!



## T-Blocks



Slimmest Possible Profile!



**5-Year Warranty!**






7625 West New York Street, Indianapolis IN 46214-4911  
 Phone: 317-273-6960 Fax: 317-273-6979  
 e-mail: [orderdesk@floteco2.com](mailto:orderdesk@floteco2.com)

[www.floteco2.com](http://www.floteco2.com)  
 Order Desk: 800-401-1723








**Y-Block  
Style**



**T-Block  
Style**



**Green** is Oxygen and **Yellow** for Medical Air, All other Medical Gases available on request

Part Number	reference	Inlet Type	Part Number	reference
T10-B70P6P6 T10-B7EP6P6 (w/ DCV)		DISS Nut & Nipple 	T10-270P6P6 T10-27EP6P6 (w/ DCV)	
TA04- B804P6P6 TA04- B8G4P6P6 (w/ DCV)			TA04- 2804P6P6 TA04- 28G4P6P6 (w/ DCV)	
T10-BC0P6P6 T10-BCEP6P6 (w/ DCV)		DISS Female Hand-Tight 	T10-2C0P6P6 T10-2CEP6P6 (w/ DCV)	
TA04- BL04P6P6 TA04- BLG4P6P6 (w/ DCV)			TA04- 2L04P6P6 TA04- 2LG4P6P6 (w/ DCV)	
T10-BG0P6P6 T10-BGEP6P6 (w/ DCV)		Chemetron 	T10-2G0P6P6 T10-2GEP6P6 (w/ DCV)	
TA04- BF04P6P6 TA04- BFG4P6P6 (w/ DCV)			TA04- 2F04P6P6 TA04- 2FG4P6P6 (w/ DCV)	
T10-BE0P6P6 T10-BEEP6P6 (w/ DCV)		Ohio / Ohmeda 	T10-2E0P6P6 T10-2EEP6P6 (w/ DCV)	
TA04- BD04P6P6 TA04- BDG4P6P6 (w/ DCV)			TA04- 2D04P6P6 TA04-2DG4P6P6 (w/ DCV)	
T10-B10P6P6 T10-B1EP6P6 (w/ DCV)		Shrade 	T10-210P6P6 T10-21EP6P6 (w/ DCV)	
TA04- B404P6P6 TA04- B4G4P6P6 (w/ DCV)			TA04- 2404P6P6 TA04-24G4P6P6 (w/ DCV)	

21 CFR 820, 803, 806, 807 - Subparts A-D, & (821 where applicable)  
 CE Mark 0843 CGA G4.1  
 IEC 62366-1:2015 ASTM-G175-03  
 European Council Directive 93/42/ECC

ISO 13485:2016  
 TGR 2002  
 MRI Conditional (selected P/N's)  
 ISO 10524:2006 (selected P/N's)

MDSAP Audited  
 BS EN 14971:2012  
 RoHS 3 Directive (EU) 2017/2012