



# Double Flowmeter Assemblies

## Y-Blocks




**All-Metal-Construction  
No Plastic!**



## T-Blocks



**Slimmest Possible Profile!**



**5-Year Warranty!**



Dual - Flowmeter from a single wall outlet.


Optional Pediatric and High Flowmeter

Optional removable Resident Nipple Adapter for CPAP of Humidifier applications

Optional DISS-1240 Check Valve, Female Ohio or Chemetron


Optional T or Y styles.

Available flows from 1/32nd of ...










**Y-Block  
Style**



**T-Block  
Style**



**Green** is Oxygen and **Yellow** for Medical Air, All other Medical Gases available on request

Part Number		Inlet Type	Part Number	
T10-B70P6P6 T10-B7EP6P6 (w/ DCV)		DISS Nut & Nipple 	T10-270P6P6 T10-27EP6P6 (w/ DCV)	
TA04- B804P6P6 TA04- B8G4P6P6 (w/ DCV)			TA04- 2804P6P6 TA04- 28G4P6P6 (w/ DCV)	
T10-BC0P6P6 T10-BCEP6P6 (w/ DCV)		DISS Female Hand-Tight 	T10-2C0P6P6 T10-2CEP6P6 (w/ DCV)	
TA04- BL04P6P6 TA04- BLG4P6P6 (w/ DCV)			TA04- 2L04P6P6 TA04- 2LG4P6P6 (w/ DCV)	
T10-BG0P6P6 T10-BGEP6P6 (w/ DCV)		Chemetron 	T10-2G0P6P6 T10-2GEP6P6 (w/ DCV)	
TA04- BF04P6P6 TA04- BFG4P6P6 (w/ DCV)			TA04- 2F04P6P6 TA04- 2FG4P6P6 (w/ DCV)	
T10-BE0P6P6 T10-BEEP6P6 (w/ DCV)		Ohio / Ohmeda 	T10-2E0P6P6 T10-2EEP6P6 (w/ DCV)	
TA04- BD04P6P6 TA04- BDG4P6P6 (w/ DCV)			TA04- 2D04P6P6 TA04-2DG4P6P6 (w/ DCV)	
T10-B10P6P6 T10-B1EP6P6 (w/ DCV)		Shrade 	T10-210P6P6 T10-21EP6P6 (w/ DCV)	
TA04- B404P6P6 TA04- B4G4P6P6 (w/ DCV)			TA04- 2404P6P6 TA04-24G4P6P6 (w/ DCV)	