

- Fitting end
- 12 - 1/8" MNPT
- 14 - 1/4" MNPT
- 22 - 1/8" FNPT
- 2S - 1/8" FNPT Swivel
- 90-HB-4 - 90 Elbow 1/4" Hose Barb
- DCV - DISS Check Valve
- DF - DISS Fitting
- HB-4 - 1/4" Hose Barb
- HB-5 - 5/16" Hose Barb
- HN - Hex Nut Female
- KN - Knurled Nut Female
- WO - Wall Outlet w/ Face Plate and 1/2" OD tube rotating pigtail

----- QD End

- Quick Disconnect Female (BS) British Standard Oxygen - B9
- Quick Disconnect Female Chemetron Carbon Dioxide - U2
- Quick Disconnect Female Chemetron Medical Air - D2
- Quick Disconnect Female Chemetron Nitrous Oxide - N2
- Quick Disconnect Female Chemetron Oxygen - C2
- Quick Disconnect Female Chemetron Vacuum (Suction) - V2
- Quick Disconnect Female Chemetron WAGD (EVAC) - E2
- Quick Disconnect Female Hansen Large Brass Oxygen - C6
- Quick Disconnect Female Hansen Large Stainless Oxygen - C4
- Quick Disconnect Female Hansen Small Brass Oxygen - H4
- Quick Disconnect Female Ohio Carbon Dioxide - Y2
- Quick Disconnect Female Ohio Medical Air - A2
- Quick Disconnect Female Ohio Nitrous Oxide - H2
- Quick Disconnect Female Ohio Oxygen - O2
- Quick Disconnect Female Ohio Vacuum (Suction) - J2
- Quick Disconnect Female Oxequip Oxygen - Z4
- Quick Disconnect Female Oxequip Medical Air - Z8
- Quick Disconnect Female Puritan Bennet Nitrous Oxide - K2
- Quick Disconnect Female Puritan Bennet Vacuum (Suction) - M2
- Quick Disconnect Female Puritan Bennet WAGD (EVAC) - B2
- Quick Disconnect Female Puritan Bennett Medical Air - R2
- Quick Disconnect Female Puritan Bennett Oxygen - P2
- Quick Disconnect Female Schrader Medical Air - T2
- Quick Disconnect Female Schrader Nitrogen - Q2
- Quick Disconnect Female Schrader Nitrous Oxide - G2
- Quick Disconnect Female Schrader Oxygen - S2
- Quick Disconnect Female Schrader Vacuum (Suction) - W2
- Quick Disconnect Female (AFNOR) French Knurled Nut Medical Air - H8
- Quick Disconnect Female (AFNOR) French Knurled Nut Nitrous Oxide - H7
- Quick Disconnect Female (AFNOR) French Knurled Nut Oxygen - H6
- Quick Disconnect Female (AFNOR) French Knurled Nut Vacuum (Suction) - H5
- Quick Disconnect Female (NIST) European Hex Nut Medical Air - D8
- Quick Disconnect Female (NIST) European Hex Nut Nitrous Oxide - D7
- Quick Disconnect Female (NIST) European Hex Nut Oxygen - D6
- Quick Disconnect Female (NIST) European Hex Nut Vacuum (Suction) - D5
- Quick Disconnect Female (SIS) Australian Knurled Nut Medical Air - K8
- Quick Disconnect Female (SIS) Australian Knurled Nut Nitrous Oxide - K7
- Quick Disconnect Female (SIS) Australian Knurled Nut Oxygen - K6
- Quick Disconnect Female (SIS) Australian Knurled Nut Vacuum (Suction) - K5
- Quick Disconnect Socket (SIS) Australian Medical Air - E8
- Quick Disconnect Socket (SIS) Australian Nitrous Oxide - E7
- Quick Disconnect Socket (SIS) Australian Oxygen - E6
- Quick Disconnect Socket (SIS) Australian Knurled (Suction) - E5
- Quick Disconnect Male Chemetron Carbon Dioxide - U1
- Quick Disconnect Male Chemetron Medical Air - D1
- Quick Disconnect Male Chemetron Nitrous Oxide - N1
- Quick Disconnect Male Chemetron Oxygen - C1
- Quick Disconnect Male Chemetron Vacuum (Suction) - V1
- Quick Disconnect Male Chemetron WAGD (EVAC) - E1
- Quick Disconnect Male Medstar Medical Air x 1/8" MNPT - L3
- Quick Disconnect Male Medstar Nitrous Oxide - L1
- Quick Disconnect Male Medstar Oxygen x 1/8" MNPT - L5
- Quick Disconnect Male Medstar Vacuum (Suction) x 1/8" MNPT - L7

- Quick Disconnect Male non-swivel Schrader Medical Air - T1
- Quick Disconnect Male non-swivel Schrader Nitrogen - Q1
- Quick Disconnect Male non-swivel Schrader Nitrous Oxide - G1
- Quick Disconnect Male non-swivel Schrader Oxygen - S1
- Quick Disconnect Male non-swivel Schrader Vacuum (Suction) - W1
- Quick Disconnect Male Ohio Carbon Dioxide - Y1
- Quick Disconnect Male Ohio Medical Air - A1
- Quick Disconnect Male Ohio Medical Air (ISO Black/White) - A1BW
- Quick Disconnect Male Ohio Nitrous Oxide - H1
- Quick Disconnect Male Ohio Oxygen - O1
- Quick Disconnect Male Ohio Oxygen (ISO White) - O1W
- Quick Disconnect Male Ohio Vacuum (Suction) - J1
- Quick Disconnect Male Ohio Vacuum (Suction) (ISO Yellow) - J1Y
- Quick Disconnect Male Ohio WAGD (EVAC) - F1
- Quick Disconnect Male Ohio WAGD (EVAC) - F2
- Quick Disconnect Male Oxequip Medical Air - Z3
- Quick Disconnect Male Oxequip Nitrous Oxide - Z1
- Quick Disconnect Male Oxequip Oxygen - Z5
- Quick Disconnect Male Oxequip Vacuum (Suction) - Z7
- Quick Disconnect Male Hansen Large Stainless Oxygen - C3
- Quick Disconnect Male Hansen Large Brass Oxygen - C5
- Quick Disconnect Male Hansen Small Brass Oxygen - H3
- Quick Disconnect Male Puritan Bennet Nitrous Oxide K1
- Quick Disconnect Male Puritan Bennet Vacuum (Suction) - M1
- Quick Disconnect Male Puritan Bennet Vacuum (Suction) (ISO Yellow) - M1Y
- Quick Disconnect Male Puritan Bennet WAGD (EVAC) - B1
- Quick Disconnect Male Puritan Bennett Medical Air - R1
- Quick Disconnect Male Puritan Bennett Medical Air (ISO Yellow) - R1Y
- Quick Disconnect Male Puritan Bennett Oxygen - P1
- Quick Disconnect Male Puritan Bennett Oxygen (ISO White) - P1W
- Quick Disconnect Male swivel Schrader Medical Air - T3
- Quick Disconnect Male swivel Schrader Nitrogen - Q3
- Quick Disconnect Male swivel Schrader Nitrous Oxide - G3
- Quick Disconnect Male swivel Schrader Oxygen - S3
- Quick Disconnect Male swivel Schrader Vacuum (Suction) - W3
- Quick Disconnect Male (BS) British Standard Medical Air - B8
- Quick Disconnect Male (BS) British Standard Nitrous Oxide - B7
- Quick Disconnect Male (BS) British Standard Oxygen - B6
- Quick Disconnect Male (BS) British Standard Vacuum (Suction) - B5
- Quick Disconnect Male (DIN) Germany Carbon Dioxide - F9
- Quick Disconnect Male (DIN) Germany Medical Air - F8
- Quick Disconnect Male (DIN) Germany Nitrogen - F4
- Quick Disconnect Male (DIN) Germany Nitrous Oxide - F7
- Quick Disconnect Male (DIN) Germany Oxygen - F6
- Quick Disconnect Male (DIN) Germany Vacuum (Suction) - F5
- Quick Disconnect Male (NIST) European Hex Nut Vacuum (Suction) - A5
- Quick Disconnect Male (NIST) European Hex Nut Oxygen - A6
- Quick Disconnect Male (NIST) European Hex Nut Nitrous Oxide - A7
- Quick Disconnect Male (NIST) European Hex Nut Medical Air - A8
- Quick Disconnect Male (BS) British Std w/Knurled Nut Vacuum (Suction) G5

----- PRODUCT GROUP:
 QD - Quick Disconnect Fitting

A1 - Quick Disconnect Male Ohio Medical Air	M1 - Quick Disconnect Male Puritan Bennet Vacuum (Suction)
A1BW - Quick Disconnect Male Ohio Medical Air (ISO Black/White)	M1Y - Quick Disconnect Male Puritan Bennet Vacuum (Suction) (ISO Yellow)
A2 - Quick Disconnect Female Ohio Medical Air	M2 - Quick Disconnect Female Puritan Bennet Vacuum (Suction)
A5 - Quick Disconnect Male (NIST) European Hex Nut Vacuum (Suction)	N1 - Quick Disconnect Male Chemetron Nitrous Oxide
A6 - Quick Disconnect Male (NIST) European Hex Nut Oxygen	N2 - Quick Disconnect Female Chemetron Nitrous Oxide
A7 - Quick Disconnect Male (NIST) European Hex Nut Nitrous Oxide	O1 - Quick Disconnect Male Ohio Oxygen
A8 - Quick Disconnect Male (NIST) European Hex Nut Medical Air	O1W - Quick Disconnect Male Ohio Oxygen (ISO White)
B1 - Quick Disconnect Male Puritan Bennet WAGD (EVAC)	O2 - Quick Disconnect Female Ohio Oxygen
B2 - Quick Disconnect Female Puritan Bennet WAGD (EVAC)	P1 - Quick Disconnect Male Puritan Bennett Oxygen
B5 - Quick Disconnect Male (BS) British Standard Vacuum (Suction)	P1W - Quick Disconnect Male Puritan Bennett Oxygen (ISO White)
B6 - Quick Disconnect Male (BS) British Standard Oxygen	P2 - Quick Disconnect Female Puritan Bennett Oxygen
B7 - Quick Disconnect Male (BS) British Standard Nitrous Oxide	Q1 - Quick Disconnect Male non-swivel Schrader Nitrogen
B8 - Quick Disconnect Male (BS) British Standard Medical Air	Q2 - Quick Disconnect Female Schrader Nitrogen
B9 - Quick Disconnect Female (BS) British Standard Oxygen	Q3 - Quick Disconnect Male swivel Schrader Nitrogen
C1 - Quick Disconnect Male Chemetron Oxygen	R1 - Quick Disconnect Male Puritan Bennett Medical Air
C2 - Quick Disconnect Female Chemetron Oxygen	R1Y - Quick Disconnect Male Puritan Bennett Medical Air (ISO Yellow)
C3 - Quick Disconnect Male Hansen Large Stainless Oxygen	R2 - Quick Disconnect Female Puritan Bennett Medical Air
C4 - Quick Disconnect Female Hansen Large Stainless Oxygen	S1 - Quick Disconnect Male non-swivel Schrader Oxygen
C5 - Quick Disconnect Male Hansen Large Brass Oxygen	S2 - Quick Disconnect Female Schrader Oxygen
C6 - Quick Disconnect Female Hansen Large Brass Oxygen	S3 - Quick Disconnect Male swivel Schrader Oxygen
D1 - Quick Disconnect Male Chemetron Medical Air	S5 - Quick Disconnect Male Scandinavian (AGA) Vacuum (Suction)
D2 - Quick Disconnect Female Chemetron Medical Air	S6 - Quick Disconnect Male Scandinavian (AGA) Oxygen
D5 - Quick Disconnect Female (NIST) European Hex Nut Vacuum (Suction)	S7 - Quick Disconnect Male Scandinavian (AGA) Nitrous Oxide
D6 - Quick Disconnect Female (NIST) European Hex Nut Oxygen	S8 - Quick Disconnect Male Scandinavian (AGA) Medical Air
D7 - Quick Disconnect Female (NIST) European Hex Nut Nitrous Oxide	T1 - Quick Disconnect Male non-swivel Schrader Medical Air
D8 - Quick Disconnect Female (NIST) European Hex Nut Medical Air	T2 - Quick Disconnect Female Schrader Medical Air
E1 - Quick Disconnect Male Chemetron WAGD (EVAC)	T3 - Quick Disconnect Male swivel Schrader Medical Air
E2 - Quick Disconnect Female Chemetron WAGD (EVAC)	T5 - Quick Disconnect Male Japanese (JIS) Vacuum (Suction)
E5 - Quick Disconnect Socket (SIS) Australian Vacuum (Suction)	T6 - Quick Disconnect Male Japanese (JIS) Oxygen
E6 - Quick Disconnect Socket (SIS) Australian Oxygen	T7 - Quick Disconnect Male Japanese (JIS) Nitrous Oxide
E7 - Quick Disconnect Socket (SIS) Australian Nitrous Oxide	T8 - Quick Disconnect Male Japanese (JIS) Medical Air
E8 - Quick Disconnect Socket (SIS) Australian Medical Air	U1 - Quick Disconnect Male Chemetron Carbon Dioxide
F1 - Quick Disconnect Male Ohio WAGD (EVAC)	U2 - Quick Disconnect Female Chemetron Carbon Dioxide
F2 - Quick Disconnect Male Ohio WAGD (EVAC)	V1 - Quick Disconnect Male Chemetron Vacuum (Suction)
F4 - Quick Disconnect Male (DIN) Germany Nitrogen	V2 - Quick Disconnect Female Chemetron Vacuum (Suction)
F5 - Quick Disconnect Male (DIN) Germany Vacuum (Suction)	W1 - Quick Disconnect Male non-swivel Schrader Vacuum (Suction)
F6 - Quick Disconnect Male (DIN) Germany Oxygen	W2 - Quick Disconnect Female Schrader Vacuum (Suction)
F7 - Quick Disconnect Male (DIN) Germany Nitrous Oxide	W3 - Quick Disconnect Male swivel Schrader Vacuum (Suction)
F8 - Quick Disconnect Male (DIN) Germany Medical Air	Y1 - Quick Disconnect Male Ohio Carbon Dioxide
F9 - Quick Disconnect Male (DIN) Germany Carbon Dioxide	Y2 - Quick Disconnect Female Ohio Carbon Dioxide
G1 - Quick Disconnect Male non-swivel Schrader Nitrous Oxide	Z1 - Quick Disconnect Male Oxequip Nitrous Oxide
G2 - Quick Disconnect Female Schrader Nitrous Oxide	Z3 - Quick Disconnect Male Oxequip Medical Air
G3 - Quick Disconnect Male swivel Schrader Nitrous Oxide	Z4 - Quick Disconnect Female Oxequip Oxygen
G6 - Quick Disconnect Male (BS) British Standard w/Knurled Nut Medical Air	Z5 - Quick Disconnect Male Oxequip Oxygen
G7 - Quick Disconnect Male (BS) British Standard w/Knurled Nut Oxygen	Z7 - Quick Disconnect Male Oxequip Vacuum (Suction)
G8 - Quick Disconnect Male (BS) British Standard w/Knurled Nut Nitrous Oxide	Z8 - Quick Disconnect Female Oxequip Medical Air
G9 - Quick Disconnect Male (BS) British Standard w/Knurled Nut Vacuum (Suction)	
H1 - Quick Disconnect Male Ohio Nitrous Oxide	
H2 - Quick Disconnect Female Ohio Nitrous Oxide	N4 - Quick Disconnect Female (DIN) Germany w/Knurled Nut Medical Air
H3 - Quick Disconnect Male Hansen Small Brass Oxygen	N5 - Quick Disconnect Female (DIN) Germany w/Knurled Nut Oxygen
H4 - Quick Disconnect Female Hansen Small Brass Oxygen	
H5 - Quick Disconnect Female (AFNOR) French Knurled Nut Vacuum (Suction)	
H6 - Quick Disconnect Female (AFNOR) French Knurled Nut Oxygen	
H7 - Quick Disconnect Female (AFNOR) French Knurled Nut Nitrous Oxide	
H8 - Quick Disconnect Female (AFNOR) French Knurled Nut Medical Air	
J1 - Quick Disconnect Male Ohio Vacuum (Suction)	
J1Y - Quick Disconnect Male Ohio Vacuum (Suction) (ISO Yellow)	
J2 - Quick Disconnect Female Ohio Vacuum (Suction)	
J6 - Quick Disconnect Male (DIN) Germany w/Knurled Nut Medical Air	
J7 - Quick Disconnect Male (DIN) Germany w/Knurled Nut Oxygen	
J8 - Quick Disconnect Male (DIN) Germany w/Knurled Nut Nitrous Oxide	
J9 - Quick Disconnect Male (DIN) Germany w/Knurled Nut Vacuum (Suction)	
K1 - Quick Disconnect Male Puritan Bennet Nitrous Oxide	
K2 - Quick Disconnect Female Puritan Bennet Nitrous Oxide	
K5 - Quick Disconnect Female (SIS) Australian Knurled Nut Vacuum (Suction)	
K6 - Quick Disconnect Female (SIS) Australian Knurled Nut Oxygen	
K7 - Quick Disconnect Female (SIS) Australian Knurled Nut Nitrous Oxide	
K8 - Quick Disconnect Female (SIS) Australian Knurled Nut Medical Air	
L1 - Quick Disconnect Male Medstar Nitrous Oxide	
L3 - Quick Disconnect Male Medstar Medical Air x 1/8" MNPT	
L5 - Quick Disconnect Male Medstar Oxygen x 1/8" MNPT	
L7 - Quick Disconnect Male Medstar Vacuum (Suction) x 1/8" MNPT	
N4 - Quick Disconnect Female (DIN) Germany w/Knurled Nut Medical Air	
N5 - Quick Disconnect Female (DIN) Germany w/Knurled Nut Oxygen	